



**Model: ES-PC1080HD**



**Commodity: Full HD endoscope video camera system**

MOQ: 1PCS

Warranty: 2-YEAR

Descriptions:

Active-1920 (H) ×1080(V)

H resolution-1080P video output

minimum Illumination-20Lx F1.4

Required Illumination-2000LxF5.6

S/N ratio-54dB

Video output-model HDMI, (1080/50p) RGB/YPbPr,S-VIDEO

**Model: ES-VICAM 2600HD**



**Commodity: Full HD integrated endoscope imaging system**

MOQ: 1PCS

Warranty: 2-YEAR

Descriptions:

26" Full HD with LED Screen

Integrated 4-in-1 endoscope imaging systems

Integrated Medical display\*1

EndoCamera\*1

Cold Light Source\*1

WorkStation\*1

**Model: ES-VICAM 2200HD**



**Commodity: Full HD integrated endoscope imaging system**

MOQ: 1PCS

Warranty: 2-YEAR

Descriptions:

22" Full HD with LED Screen

Integrated 4-in-1 endoscope imaging systems

Integrated Medical display\*1

EndoCamera\*1

Cold Light Source\*1

WorkStation\*1

**Model: ES-VICAM 1200HD**



**Commodity: Full HD integrated endoscope imaging system**

MOQ: 1PCS

Warranty: 2-YEAR

Descriptions:

- 12" Full HD with LED Screen
- Integrated 4-in-1 endoscope imaging systems
- Integrated Medical display\*1
- EndoCamera\*1
- Cold Light Source\*1
- WorkStation\*1

**Model: ES-VI168HD**



**Commodity: Full HD endoscope camera system**

MOQ: 1PCS

Warranty: 2-YEAR

Descriptions:

- ES-Vi168HD is self-owned researched and developed high quality endoscope cameras system.
- Adopted original SONY 1/2.8 CMOS with 3840\*2160 high definition that reached 1200lines outlet images.
- Real-time video 60 FPS that made the camera system high performance.

**Model: ES-CL300**



**Commodity: CCD camera integrated light source system**

MOQ: 1PCS

Warranty: 1-YEAR

Descriptions:

- Portable CCD Camera integrated with Light Source
- Adopt high-end CCD and LED as the Video & Light Source systems;
- Simple, easy to carry;
- Specialized and used in medical endoscopic animation, picture display and collection

**Model: ES-LC250**



**Commodity: Xenon cold light source**

MOQ: 1PCS

Warranty: 1-YEAR

Descriptions:

- Compatible with WOLF, STORZ, STRYKER & ACMI Etc.;
- Color Temperation: 7,000K
- 6 Illumination: 5.8×10 IX
- Input Voltage: AC 90~240v
- Lamp Power: 250W
- Lamp Life: 500h

**Model: ES-LC200**



**Commodity: 80W LED cold light source**

MOQ: 1PCS

Warranty: 1-YEAR

Descriptions:

Supply voltage: 100~240V 50Hz;

Power $\leq$ 250VA;

Lamp rated power  $\leq$ 80VA;

Luminous flux: 100m;

Lamp life  $\geq$ 20,000h;

Color temp: 3,000K~7,000K;

Illuminance  $\geq$  1,800,000LX;

Ra $\geq$ 90

**Model: ES-LC300**



**Commodity: 80W LED cold light source**

MOQ: 1PCS

Warranty: 1-YEAR

Descriptions:

Supply voltage: 100~240V 50Hz;

Power $\leq$ 250VA;

Lamp rated power  $\leq$ 80VA;

Luminous flux: 100m;

Lamp life  $\geq$ 20,000h;

Color temp: 3,000K~7,000K;

Illuminance  $\geq$  1,800,000LX;

Ra $\geq$ 90

**Model: ES-CI800**



**Commodity: CO2 insufflator**

MOQ: 1PCS

Warranty: 2-YEAR

Descriptions:

Supply voltage: 110~240V 50Hz

Power: 40VA

Range of flow: 2L/min~40L/min

Range of pressure: 5~25mmHg

Noise $\leq$ 50dB(A)

**Model: ES-IP700**



**Commodity: Integrated irrigation and suction system**

MOQ: 1PCS

Warranty: 2-YEAR

Descriptions:

Supply Voltage: 110~240V 50Hz

Power: 150VA

Pressure range of irrigation: 50~700mmHg

Flow range of irrigation: 10~1,500mL/min

Pressure range of suction: 15~400mmHg

Flow range of suction: 1,000~30,000mL/min

Noise $\leq$ 70dB(A)

**Model: ES-TS616**



**Commodity: Trolley of Endoscope Systems**

MOQ: 1PCS

Warranty: 2-YEAR

Descriptions:

1. Color: Shining silver
2. Material: Aluminum Alloy and Steel
3. With adjustable shelf, shelf numbers could be adjusted according to the convenience;
4. With 2 endoscopes holding device, and 1 coupler holding device;
5. With flexible monitor mount;
6. With flexible shelf for holding the computer key board and mouse;
7. Easy assembly and shipment;
8. 2 packages of cartons box,
9. Dimensions are: 710\*700\*285mm 1145\*305\*235mm
10. GW: 50KG

**Model: ESB-05**



**Commodity: MONO/BIPOLAR ESU with ARGON option**

MOQ: 1PCS

Warranty: 1-YEAR

Descriptions:

Supply voltage: 110~240V 50Hz Power: 40VA

MAX POWER 400 WATTS

7 MONOPOLAR OUTPUT AND 3 BIPOLAR OUTPUT

9 PROGRAMMABLE MEMORY

AUTOMATIC ALIMENTATION SOURCE COMPENSATION AND TISSUE VARIATION CONTROL

3 DIFFERENT INSTRUMENTS CAN BE SIMULTANEOUSLY CONFIGURED

OPTIONAL ADAPTER FOR LAPAROSCOPIC AND THORACOSCOPY

APPLICATION

ARGON MODULE OPTIONAL

TROLLEY OPTIONAL

**Model: ESS-MAXGEN**



**Commodity: MONO/BIPOLAR ESU**

MOQ: 1PCS

Warranty: 1-YEAR

Descriptions:

LCD color touchscreen

Max power 400 W

100 memory programs

Endoscopic cut modes



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail: [danielec@esse3-medical.com](mailto:danielec@esse3-medical.com)  
web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

