



25591

Digital Thermometer



LITTLE ACCESSORIES



25591

Description

Multifunctional infrared non contact thermometer for body and objects.

Effective distance 5-15 cm.

Body temperature with fever teller.

Wide range from 32.0 and 42.9 °C and green, yellow, red backlight.

Object temperature 0-100 °C (32-212 °F)

Body object function for ambient temperature, milk or bath surface temperature.

- Fever alarm
- 32 memories
- Instant response time: 1 second
- Auto/manual power off.
- Supplied with 2AA batteries.
- Size 160x95x46 mm.
- Display 30x31 mm.
- Reading scale °C/°F switchable
- Packed in multilanguage carton
- Multilanguage box and manual: GB, FR, IT, ES.



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

