



DATA SHEET

Model: 344xx Line

HEINE F.O. Laryngoscopes




MILLER



PAED (paediatric)



MC-INTOSH



HEINE CLASSIC+[®] F.O. LARYNGOSCOPES

WE introduce the "+" range of F.O. blades with HEINE HiLite[®] Fiber Optic bundle with over 6,000 fibre optic and the state of the art HEINE XHL[®] Xenon Halogen bulbs or LED HQ in the Handle. The innovative HEINE HiLite[®] Fiber Optic bundle with up to 100% more cross-section improves light transmission and makes the blade more resistant to autoclaving than ever. HEINE offers a 5-year guarantee or 4000 standard autoclave cycles of 134°C/5 min. with more than 1000 lux on all blades with HEINE HiLite[®] fibers. HEINE XHL[®] Xenon Halogen bulbs incorporate the latest Xenon Halogen Technology to give you up to 40% more light! All HEINE F.O.+ laryngoscope systems offer exceptional value for money and effectively save money for you. HEINE XHL[®] Xenon Halogen 3.5 V bulbs last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscopes bulbs. HEINE 3.5 V systems use a long-lasting 3.5 V rechargeable battery which avoids the need to replace non-rechargeable batteries at frequent intervals. HEINE F.O. laryngoscope blades and handles are compatible with the ISO 7376 (green standard).

Model	HEINE FIBRE OPTIC BLADES	Overall length	Distal length
Cod. 34411	Miller n° 00	76 mm	9.5 mm
Cod. 34412	Miller n° 0	80 mm	10.5 mm
Cod. 34413	Miller n° 1	100 mm	11.5 mm
Cod. 34414	Miller n° 2	152 mm	13 mm
Cod. 34415	PAED 0	82 mm	12 mm
Cod. 34416	PAED 1	95 mm	13 mm
Cod. 34417	Mc-Intosh n° 0	82 mm	9 mm
Cod. 34418	Mc-Intosh n° 1	93 mm	11 mm
Cod. 34419	Mc-Intosh n° 2	115 mm	13 mm
.... Cod. 34420	Mc-Intosh n° 3	135 mm	15 mm
Cod. 34421	Mc-Intosh n° 4	155 mm	15 mm
Cod. 34422	Mc-Intosh n° 5	176 mm	15 mm
Cod. 34423	MAC 3 - flex tip+	135 mm	15 mm
Model	HEINE LARYNGOSCOPE HANDLES AND CHARGER		
Cod. 34400	Short F.O. handle 2.5V for pediatric and obese		
Cod. 34401	Standard F.O. handle 2.5V		
Cod. 34405	Standard F.O. handle 2.5 LED		
Cod. 34402	Standard handle 3.5V F.O. LED Li-ion L – rechargeable - 6,500 Kelvin - 12,000 LUX at 40 mm - 50,000 hours bulb life - Handle with charge, status display.		
Cod. 31769	Spare bulb for 34400 and 34401		
Cod. 34448	Re-chargeable battery for 34402		
Cod. 34449	NT 300 charger unit for 2 handles 14x6x6.6 cm Two independent charging compartments, Fast charge, Automatic charging, Charger indicator		
Model	HEINE CLASSIC+ [®] FIBRE OPTIC LARYNGOSCOPE SETS*		
Cod. 34403	3 blades set with paed 1, mac 2, mac 3		
Cod. 34404	4 blades set with paed 1, mac 2, mac 3, mac 4		
* set includes spare bulb and F.O. handle 2.5 V (34401)			



 34411	 34412	 34413	 34414
 34415	 34416	 34417	 34418
 34419	 344120	 344121	 344122
 344123	 31769	 34448	 34449
 34400	 34401	 34405	 34402
 34403	 34404		