



# OTOSCOPE



**Cod. ESF8-680200**



## OTOSCOPE

Compact pocket otoscope Direct illumination. Enhanced Xenon Halogen Technology with 100% more light compared with conventional bulbs. Easy to maintain, flexible. Viewing window with 3x magnification and optimized casing surface for razor-sharp images and minimal reflection. Swivelling viewing window built into instrument. Useful for instrumentation, cannot be mislaid. Attachment clip with integrated on/off switch.



## OTOSCOPE

Otoscope de poche compact. Illumination directe. Technologie Xénon Halogène, qui apporte 100% de lumière en plus comparée avec les ampoules classiques. Utilisation simple, Fenêtre d'observation avec un grossissement de 3 fois et surface du boîtier optimisée pour des images des plus précises et des reflets minimisés. Fenêtre d'observation intégrée mais pouvant pivoter sur le côté. Optimal pour intervention avec grossissement. Clip de fixation avec interrupteur intégré. Bonne prise en main de l'instrument. Résistant aux chocs, robuste, anti-dérapant



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

