



ESD-369D

Micro Surgical Drill and
Saw serie



OPERATING THEATRE

**ESD-369D.500****Micro Cannulated Drill**

For Hand Surgery

Standard configuration:

- One hand piece
- Two battery 7.2V
- One charger
- One drill chuck and 1 key
- Two aseptic battery transfer rings
- One aluminum box

Specifications

- Non-load noise ≤ 65 dB
- Output power 95W
- Torque 6.5Nm
- Non-load speed 0-1000rpm

**ESD-369D.501****Micro Oscillating Saw**

For Hand Surgery

Standard configuration:

- One hand piece
- Two battery 7.2V
- One charger
- Two saw blades
- Two aseptic battery transfer rings
- One aluminum box

Specifications

- Non-load speed 0-18000osc/min
- Non-load noise ≤ 45 dB (A)
- Diameter run-out ≤ 0.5 mm
- Output Power 50W

**ESD-369D.502****Micro Bone Drill**

For Hand Surgery (Non-Cannulated)

Standard configuration:

- One hand piece
- Two battery 7.2V
- One charger
- One drill chuck and 1 key
- Two aseptic battery transfer rings
- One aluminum box

Specifications

- Non-load speed 0-1100r/min
- Non-load noise 40dB
- Diameter run-out ≤ 0.5 mm
- Output Power 50W



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

