




DATA SHEET 

Model: EST-A100

Electric Operating Table



- Electric operating table is made specially to meet demand of general surgical operations. It adopts full 304 stainless steel, excellent properties on anti-rust.
- Tabletop is X-raying available. can be sliding for over 300mm, that meet C-arm need.
- Length/Width: 2070mm/550mm
- Elevation (up/down): 700mm/1000mm
- Trendelenburg: $\geq 25^\circ$
- Reverse Trendelenburg: $\geq 25^\circ$
- Lateral Tilt (left/right): $25^\circ/25^\circ$
- Head Rest (up/down): $45^\circ/90^\circ$
- Footplate (up/down & horizontal positioning): 90°
- Back plate (up/down): $75^\circ/20^\circ$
- Kidney bridge: 120mm
- Sliding: 300mm
- Power: 220v, 50Hz, 1.0kW

Standard accessories:

- | | |
|-----------------------|------------------|
| - 2 Shoulder Support | - 2 Foot plate |
| - 2 Arm Support | - 8 Fixing Clamp |
| - 2 Body Support | - 1 Mattress |
| - 2 Long Fixing Clamp | - 1 Hand remote |
| - 2 Knee Crutch | |



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

