



ESF5-35/xx

Minor surgery suction pump



OPERATING THEATRE



ESF5-35/xx

Minor surgery suction Pump

Portable suction machine for professional use and surgery treatment, with powerful piston motor. ABS anti-shock case, light, easy to carry and keep clean. Supplied with 1 lt. jar sterilizable, with overflow valve. Complete with suction regulator, vacuum gauge, disposable filter and set of tubes.



ESF5-35/00

Technical Features

Motor	Maintenance-free Oil free piston pump
Weight	4,5 kg
Power supply	230V/50-60Hz (0,260A)
Power consumption	60 W
Suction power	70 l/min
Max vacuum	0,80 bar - 600 mmHg / -80 kPa
Operating cycle	60min ON / 30 min OFF
Noisiness	61,7 Db

Available Configurations

Jar	230V 50-60Hz	110V 60Hz
1 x 1 L	ESF5-35/00	ESF5-35/11
1 x 2 L	ESF5-35/20	ESF5-35/31
2 x 1 L	ESF5-35/40	ESF5-35/41
2 x 2 L	ESF5-35/60	ESF5-35/71

Accessories On Request

ESF5-35/xx/P	Footswitch predisposition
ESF5-35/xx/PF	Footswitch predisposition and footswitch

ESF5-35 VERSIONS ON CASTORS

Portable suction machine for different applications with powerful piston motor. Mounted on metallic base with 4 castors, two of which with brakes.

Available Configurations

Jar	230V 50-60Hz	110V 60Hz
1 x 1 L	ESF5-35/05	ESF5-35/06
2 x 1 L	ESF5-35/45	ESF5-35/46
2 x 2 L	ESF5-35/65	ESF5-35/66
2 x 4 L	ESF5-35/95	ESF5-35/97



ESF5-35/65



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

