

Ophthalmic Equipment -Portable Slit Lamp

ESLM-6 L / 6 M / 6 H Portable slit lamp



Feature

1. Portable slit lamp adopts high-brightness, cold-light illuminator, without external power supply.
2. High-brightness LED, cold-light illuminator, its brightness up to 30000 LUX, service life more than 50000 hours.
3. Sharp slit, the width can be adjusted continuously.
4. Telescopic forehead-rest, easier to observe and adjust the diopter.
5. Easier to adjust the magnification by rotating objective lens.
6. Small and light, only 700g including battery.
8. Full power supplied by battery can drive the slit lamp for 5-6 hours.
9. Iphone or Sumsung mobile adaptor optional.can take images by your mobile.

Item / Model	ESLM-6L	ESLM-6	ESLM-6H	
Microscope	Microscope Type			
	Objective Lens	1.0X	1.0X 1.6X Converging	1.0x 1.6x
	Eyepiece	10X	10X	10x 16X
	Total Magnifications	10X	10X 16X	10x 16x 25X
	Working Distance	81mm		
	Intraocular Distance	45mm-78mm		
	Diopter Adjustment	-D8 ~ +D8		
	Illumination	Slit Width		
Slit Height		1~12mm		
Slit Angle		Horizontal+30°		
Light Spot Diameter		Φ0.2mm ;Φ1 mm; Φ5mm; Φ12mm		
Filters		Heat absorption; Yellow;Redfree (Green); Cobalt Blue		
Illumination Bulb		High-Brightness White LED		
Brightness		30000Lux		
Power		5VA		
Power	Charging Voltage			
	LED Working Voltage	DC 3V, 350mAh		
	Battery	Li-ion Battery; 7.4V, 680mAh, Full Battery can supply 5-6 hours continuously		
Weight and Dimension	Dimension			
	Weight	169.5 ×109.5 ×228.5mm 700g (Including Battery)		



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

