



Motorized tables

AT series



OPHTHALMOLOGY / ENT

AT series

Motorized tables

Models



AT-16

Table size: 480*560mm

Table rise height: 630-810mm

In put power: 220V/50Hz 110V/60Hz 160W

Weight: 26KG



AT-20

High quality ,Stability movement,low noisy and easy operation.

Table size: 56*48 cm

Table elevator height: 63-82cm.

Max loading: 50~75KG

Input: 220V/50hz 110v/60Hz

weight: 18kg



AT-110

Electric table

Table Size: 480 X 540mm

Range of elevation: 630mm - 820mm

Max Loading: 55kg

Table up and down range: 180mm

Electric Voltage: 220V/50Hz or 110V/60Hz

Power: 120W

Weight: 19kg

Volume: 60.8(L)x53(w)x26.5(H)cm